

Anti-BANK1 Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP2168-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to BANK1		
Immunogen	Recombinant protein of human BANK1		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human BANK1 protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	BANK1		
Alternative Names	B-cell scaffold protein with ankyrin repeats		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
55024			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q8NDB2			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP2168-C; Recommmended standard protein: CST2168		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com