


# Product Data Sheet

## Human NUDT21 Knockdown HEK293T Cell Lysate


| Catalog #                 | Source                                                                                                                                                                                                                                                                                                                                                                                                                                             | Reactivity | Applications |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CKD0208                   | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            | Human      |              |
| <b>Description</b>        | Human NUDT21 knockdown HEK293T cell lysate is engineered by our optimized transduction of the specific shRNA.                                                                                                                                                                                                                                                                                                                                      |            |              |
| <b>Form</b>               | Liquid in 50 mM Tris-HCl, 2% SDS, 10% glycerol, 100 mM DTT, 0.01% Bromophenol blue                                                                                                                                                                                                                                                                                                                                                                 |            |              |
| <b>Gene Symbol</b>        | NUDT21                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |              |
| <b>Alternative Names</b>  | CFIM25; CPSF25; CPSF5; Cleavage and polyadenylation specificity factor subunit 5; Cleavage and polyadenylation specificity factor 25 kDa subunit; CFIm25; CPSF 25 kDa subunit; Nucleoside diphosphate-linked moiety X motif 21; Nudix motif 21; Pre-mRNA cleava                                                                                                                                                                                    |            |              |
| <b>Entrez Gene</b>        | 11051                                                                                                                                                                                                                                                                                                                                                                                                                                              |            |              |
| <b>SwissProt</b>          | O43809                                                                                                                                                                                                                                                                                                                                                                                                                                             |            |              |
| <b>KO/KD Validation</b>   | RT-qPCR                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Cell Line</b>          | HEK293T                                                                                                                                                                                                                                                                                                                                                                                                                                            |            |              |
| <b>Components</b>         | Human NUDT21 knockdown HEK293T cell lysate x 1 vial; shControl HEK293T cell lysate x 1 vial                                                                                                                                                                                                                                                                                                                                                        |            |              |
| <b>Directions for Use</b> | Thaw the vial on ice or at room temperature. Mix gently by pipetting or brief vortexing. Briefly centrifuge to collect the sample. Load the desired amount of lysate directly onto an SDS-PAGE gel. Perform SDS-PAGE and Western blot analysis using standard laboratory protocols. Note: This lysate has been pre-denatured in SDS sample buffer containing DTT and does not require additional reducing agent or boiling before electrophoresis. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Upon delivery aliquot and store at -80°C for one year. Avoid freeze/thaw cycles.                                                                                                                                                                                                                                                                                                                                                   |            |              |

### COHESION BIOSCIENCES LIMITED


# Product Data Sheet




RT-qPCR analysis of NUDT21 knockdown efficiency. Compared with the 293T shRNA control (left), NUDT21 mRNA expression in the stable shRNA knockdown cell line (right) decreased, and the knockdown efficiency was 95%.



RT-qPCR amplification curves. The Ct values of the knockout sample significantly increased, confirming the effective knockdown of NUDT21 mRNA. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



RT-qPCR melt curve analysis demonstrated highly specific PCR amplification product without detectable nonspecific amplification or primer-dimer formation. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]



This picture presents the raw fluorescence data (relative fluorescence units, RFU) used to generate the melt curves. [Red curve (knockdown sample), Blue curve (shRNA control), Purple and Green curve (reference gene GAPDH)]

**COHESION BIOSCIENCES LIMITED**