

## Human IgA

Catalog #	Source	Reactivity	Applications
CRP9034	Recombinant Protein		WB, IH, IP

**Description** The purified total Human IgA can be used as a negative control for antibodies in Western blot, immunoprecipitation, and immunohistochemistry, etc.

**Form** Purified Human IgA in PBS, pH 7.3, 0.01% sodium azide.

**Gene Symbol**

**Alternative Names**

**Entrez Gene**

**SwissProt**

**Purity** Greater than 95% as determined by reducing SDS-PAGE.

**Chemical Structure**

**Quality Control**

**Directions for Use** Please minimize freeze-thaw cycles.

**Storage/Stability** Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**[www.cohesionbio.com](http://www.cohesionbio.com)**ORDER**[order@cohesionbio.com](mailto:order@cohesionbio.com)**SUPPORT**[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)**CUSTOM**[custom@cohesionbio.com](mailto:custom@cohesionbio.com)