

## Recombinant HPV16 L1 Protein

Catalog #	Source	Reactivity	Applications
CRP9109	Recombinant Protein	HPV	ELISA, LFIA, CLIA

**Description** Recombinant HPV16 L1 Protein is produced by E.coli expression system.

**Form** Liquid in a 0.2 µm filtered solution of 20mM Tris, pH 8.0.

**Gene Symbol**

**Alternative Names**

**Entrez Gene**

**SwissProt**

**Purity** Greater than 95% as determined by reducing SDS-PAGE.

**Chemical Structure**

**Quality Control**

**Directions for Use** Please minimize freeze-thaw cycles.

**Storage/Stability** Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

## WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

## ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

## SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

## CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)