

Recombinant Human PSME1 protein

Catalog #	Source	Reactivity	Applications
CST1081	E. coli	H	
Description	Recombinant Human PSME1 protein is produced by E. coli expression system and the target gene encoding Met1-Tyr249 is expressed with N-terminus GST tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PSME1		
Alternative Names	IFI5111; Proteasome activator complex subunit 1; 11S regulator complex subunit alpha; REG-alpha; Activator of multicatalytic protease subunit 1; Interferon gamma up-regulated I-5111 protein; IGUP I-5111; Proteasome activator 28 subunit alpha; PA28a; PA28a		
Entrez Gene	5720		
SwissProt	Q06323		
Purity			
Chemical Structure	MAMLRVQPEA QAKVDVFRED LCTKTENLLG SYFPKKISEL DAFLKEPALN EANLSNLKAP LDIPVPDPVK EKEKEERKKQ QEKEDKDEKK KGEDEDKGPP CGPVNCNEKI VLLQRLKPE IKDVIEQLNL VTTWLQLQIP RIEDGNNFGV AVQEKVFELM TSLHTKLEGF HTQISKYFSE RGDAVTAAK QPHVGDYRQL VHELDEAEYR DIRLMVMEIR NAYAVLYDII LK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com