

Recombinant Human TRPS1 Protein

Catalog #	Source	Reactivity	Applications
CST1129	E. coli	H	
Description	Recombinant Human TRPS1 Protein is produced by E. coli expression system and the target gene encoding Asn11-Gln220 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRPS1		
Alternative Names	Zinc finger transcription factor Trps1; Tricho-rhino-phalangeal syndrome type I protein; Zinc finger protein GC79		
Entrez Gene	7227		
SwissProt	Q9UHF7		
Purity			
Chemical Structure	NVASEGEGQI LEPIGTESKV SGKNKEFSAD QMSENTDQSD AAELNHKEEH SLHVQDPSSS SKKDLKSAVL SEKAGFNYES PSKGGNFPSF PHDEVTDNRNM LAFSSPAAGG VCEPLKSPQR AEADDPQDMA CTPSGDSLET KEDQKMSPKA TEETGQAQSG QANCQGLSPV SVASKNPQVP SDGGVRLNKS KTDLLVNDNP DPAPLSPELQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com