

Product Data Sheet

Recombinant Human ZFP90 Protein

Catalog #	Source	Reactivity	Applications
CST1130	E. coli	H	
Description	Recombinant Human ZFP90 Protein is produced by E. coli expression system and the target gene encoding Asp101-Glu230 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZFP90		
Alternative Names	KIAA1954; ZNF756; Zinc finger protein 90 homolog; Zfp-90; Zinc finger protein 756		
Entrez Gene	146198		
SwissProt	Q8TF47		
Purity			
Chemical Structure	DVSEVSHCTH DLLHATLEDS WDVSSQLDRQ QENWKRHLGS EASTQKKIIT PQENFEQNKFG GENSRLNTNL VTQLNIPARI RPSECETLGS NLGHNADLLN ENNILAKKKP YKCDKCRKAF IHRSSLTKHE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com