

Product Data Sheet

Recombinant Human HRPT2 Protein

Catalog #	Source	Reactivity	Applications
CST1135	E. coli	H	
Description	Recombinant Human HRPT2 Protein is produced by E. coli expression system and the target gene encoding Lys331-Asp410 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDC73		
Alternative Names	C1orf28; HRPT2; Parafibromin; Cell division cycle protein 73 homolog; Hyperparathyroidism 2 protein		
Entrez Gene	79577		
SwissProt	Q6P1J9		
Purity			
Chemical Structure	KTQTPAAQPV PRPVSQARPP PNQKKGSRTP IIIIPAATTS LITMLNAKDL LQDLKFVPSD EKKKQGCQRE NETLIQRRKD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com