

Product Data Sheet

Recombinant Human MAMLD1 Protein

Catalog #	Source	Reactivity	Applications
CST1137	E. coli	H	
Description	Recombinant Human MAMLD1 Protein is produced by E. coli expression system and the target gene encoding Ser491-Gln600 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MAMLD1		
Alternative Names	CG1; CXorf6; Mastermind-like domain-containing protein 1; F18; Protein CG1		
Entrez Gene	10046		
SwissProt	Q13495		
Purity			
Chemical Structure	SQQQQQQQQQ QQANVIFKPI SSNSSKTL SMIMQQGMASSS PGATEPFTFG NTKPLSHFVS EPGPQKMPSM PTTSRQPSLL HYLQQPTPTQ ASSATASSTA TATLQLQQQQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com