

Product Data Sheet

Recombinant Human BCL2L2 Protein

Catalog #	Source	Reactivity	Applications
CST1141	E. coli	H	
Description	Recombinant Human BCL2L2 Protein is produced by E. coli expression system and the target gene encoding Arg11-Val170 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	BCL2L2		
Alternative Names	BCLW; KIAA0271; Bcl-2-like protein 2; Bcl2-L-2; Apoptosis regulator Bcl-W		
Entrez Gene	599,100529063		
SwissProt	Q92843		
Purity			
Chemical Structure	RALVADFGY KLRQKGYVCG AGPGEGPAAD PLHQAMRAAG DEFETRFRT FSDLAAQLHV TPGSAQQRFT QVSDELFQGG PNWGRLVAFV VFGAALCAES VNKEMEPLVG QVQEWMVAYL ETQLADWIHS SGGWAEFTAL YGDGALEEAR RLREGNWASV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com