

Product Data Sheet

Recombinant Human UBE2M Protein

Catalog #	Source	Reactivity	Applications
CST1144	E. coli	H	
Description	Recombinant Human UBE2M Protein is produced by E. coli expression system and the target gene encoding Lys11-Phe80 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UBE2M		
Alternative Names	UBC12; NEDD8-conjugating enzyme Ubc12; NEDD8 carrier protein; NEDD8 protein ligase; Ubiquitin-conjugating enzyme E2 M		
Entrez Gene	9040		
SwissProt	P61081		
Purity			
Chemical Structure	KKEESAGGT KGSSKKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLNF KLVICPDEGF YKSGKFVFSF		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com