

Product Data Sheet

Recombinant Human COMMD6 Protein

Catalog #	Source	Reactivity	Applications
CST1149	E. coli	H	
Description	Recombinant Human COMMD6 Protein is produced by E. coli expression system and the target gene encoding Met1-Val85 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	COMMD6		
Alternative Names	COMM domain-containing protein 6		
Entrez Gene	170622		
SwissProt	Q7Z4G1		
Purity			
Chemical Structure	MEASSEPLD AKSDVTNQLV DFQWKLGMAY SSDTCRSLKY PYVAVMLKVA DHSGQVKTKC FEMTIPQFQN FYRQFKEIAA VIETV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com