

Product Data Sheet

Recombinant Human CLTCL1 Protein

Catalog #	Source	Reactivity	Applications
CST1154	E. coli	H	
Description	Recombinant Human CLTCL1 Protein is produced by E. coli expression system and the target gene encoding Asp1531-Glu1640 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CLTCL1		
Alternative Names	CLH22; CLTCL; CLTD; Clathrin heavy chain 2; Clathrin heavy chain on chromosome 22; CLH-22		
Entrez Gene	8218		
SwissProt	P53675		
Purity			
Chemical Structure	DHLYKDAMQH AAESRDAELA QKLLQWFLEE GKRECFAACL FTCYDLLRPD MVLELAWRHN LVDLAMPYFI QVMREYLSKV DKLDALESLR KQEEHVTEPA PLVFDFDGHE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com