

Product Data Sheet

Recombinant Human SENP1 Protein

Catalog #	Source	Reactivity	Applications
CST1167	E. coli	H	
Description	Recombinant Human SENP1 Protein is produced by E. coli expression system and the target gene encoding His291-Gln440 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SENP1		
Alternative Names	Sentrin-specific protease 1; Sentrin/SUMO-specific protease SENP1		
Entrez Gene	29843		
SwissProt	Q9P0U3		
Purity			
Chemical Structure	HPHHHHSVPH QPDNLAASNT QSEGSDSVIL LKVKDSQTPT PSSTFFQAEL WIKELTSVYD SRARERLRQI EEQKALALQL QNQRLQERH SVHDSVELHL RVPLEKEIPV TVVQETQKKG HKLTDESEDEF PEITEEMEKE IKNVFRNGNQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com