

Product Data Sheet

Recombinant Human ZNF385B Protein

Catalog #	Source	Reactivity	Applications
CST1177	E. coli	H	
Description	Recombinant Human ZNF385B Protein is produced by E. coli expression system and the target gene encoding Asn131-Glu360 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZNF385B		
Alternative Names	ZNF533; Zinc finger protein 385B; Zinc finger protein 533		
Entrez Gene	151126		
SwissProt	Q569K4		
Purity			
Chemical Structure	NTMDPVQKAV INHTFGVSIP PKKKQVISCN VCQLRFNSDS QAEAHYKGSK HAKKVKALDA TKNKPKMVPS KDSAKANPSC SITPITGNNS DKSEDKGKLK ASSSSQPSSS ESGSFLLKSG TTPLPPGAAT SPSKSTNGAP GTVVESEEEK AKKLLYCSLC KVAVNSLSQL EAHNTGSKHK TMVEARNGAG PIKSYPRPGS RLKMQNGSKG SGLQNKTFHC EICDVHVNSE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com