

Product Data Sheet

Recombinant Human NAP1L3 Protein

Catalog #	Source	Reactivity	Applications
CST1180	E. coli	H	
Description	Recombinant Human NAP1L3 Protein is produced by E. coli expression system and the target gene encoding Glu191-Gly490 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NAP1L3		
Alternative Names	BNAP; Nucleosome assembly protein 1-like 3		
Entrez Gene	4675		
SwissProt	Q99457		
Purity			
Chemical Structure	EEEENPKENP EVKAEKEVP KEIPEVKDEE KEVPKEIPEV KAEKADSKD CMEATPEVKE DPKEVPQVKA DDKEQPKATE AKARAAVRET HKRVPEERLQ DSVDLKRARK GKPKREDPKG IPDYWLIVLK NVDKLGPMIQ KYDEPILKFL SDVSLKFSKP GQPVSYTFEF HFLPNPYFRN EVLVKTYIIK AKPDHNDPFF SWGWEIEDCK GCKIDWRRGK DVTVTQTQSR TT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com