

Product Data Sheet

Recombinant Human IKZF4 Protein

Catalog #	Source	Reactivity	Applications
CST1182	E. coli	H	
Description	Recombinant Human IKZF4 Protein is produced by E. coli expression system and the target gene encoding Gly431-Val520 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	IKZF4		
Alternative Names	KIAA1782; ZNFN1A4; Zinc finger protein Eos; Ikaros family zinc finger protein 4		
Entrez Gene	64375		
SwissProt	Q9H2S9		
Purity			
Chemical Structure	GGPLLYRPRG PLTDPGASPS NGCQDSTDTE SNHEDRVAGV VSLPQGPPPQ PPPTIVVGRH SPAYAKEDPK PQEGLLRGTP GPSKEVLRVV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com