

Product Data Sheet

Recombinant Human NOP2 Protein

Catalog #	Source	Reactivity	Applications
CST1187	E. coli	H	
Description	Recombinant Human NOP2 Protein is produced by E. coli expression system and the target gene encoding Gln621-Lys800 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NOP2		
Alternative Names	NOL1; Putative ribosomal RNA methyltransferase NOP2; Nucleolar protein 1; Nucleolar protein 2 homolog; Proliferating-cell nucleolar antigen p120; Proliferation-associated nucleolar protein p120		
Entrez Gene	4839		
SwissProt	P46087		
Purity			
Chemical Structure	QPAKKAKGAA KTKQQLQKQQ HPKKASFQKL NGISKGADSE LSTVPSVTKT QASSSFQDSS QPAGKAEGIR EPKVTGKLKQ RSPKLQSSKK VAFLRQNAPP KGTDTPQTPAV LSPSKTQATL KPKDHHQPLG RAKGVEKQQL PEQPFEEKAAF QKQNDTPKGP QPPTVSPIRS SRPPPAKRKK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com