

Product Data Sheet

Recombinant Human ZC3HC1 Protein

Catalog #	Source	Reactivity	Applications
CST1189	E. coli	H	
Description	Recombinant Human ZC3HC1 Protein is produced by E. coli expression system and the target gene encoding Ala341-Met480 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZC3HC1		
Alternative Names	NIPA; Nuclear-interacting partner of ALK; Nuclear-interacting partner of anaplastic lymphoma kinase; hNIPA; Zinc finger C3HC-type protein 1		
Entrez Gene	51530		
SwissProt	Q86WB0		
Purity			
Chemical Structure	AEKSPGPIVS RTRSWDSSSP VDRPEPEAAS PTTRTRPVTR SMGTGDTPGL EVPSSPLRKA KRARLCSSSS SDTSSRSFFD PTSQHRDWCP WVNITLGKES RENGTEPDA SAPAEPGWKA VTILLAHKQ SSQPAETDSM		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com