

Product Data Sheet

Recombinant Human C6orf132 Protein

Catalog #	Source	Reactivity	Applications
CST1196	E. coli	H	
Description	Recombinant Human C6orf132 Protein is produced by E. coli expression system and the target gene encoding Glu451-Pro620 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C6orf132		
Alternative Names	Uncharacterized protein C6orf132		
Entrez Gene	647024		
SwissProt	Q5T0Z8		
Purity			
Chemical Structure	ENTASSAPVD WRDPSQMEKL RNELAAYLCG SRREDRFLSH RPGPTVAPQS KEGKKGPRLP EKETLLSLPA KDTPPGVPEK SLGGSSLTET EAAPSLTLPS VDYIPQDSPT PSVRQIRNEL EARLSSAAEK EAKPSIGSLP PKPRLEGGRI CENGADDDKL SKPVAKNLPP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com