

Product Data Sheet

Recombinant Human C18orf63 Protein

Catalog #	Source	Reactivity	Applications
CST1197	E. coli	H	
Description	Recombinant Human C18orf63 Protein is produced by E. coli expression system and the target gene encoding Glu371-Gln660 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C18orf63		
Alternative Names	Uncharacterized protein C18orf63		
Entrez Gene	644041		
SwissProt	Q68DL7		
Purity			
Chemical Structure	ELSVSQPTSG IFSALHLQPE SVQGRKKSLS IRAPQVHSEV LMPNRRGNTQV QHTNLSSQSN ITPKFVPVFK NRLLQMNKNT SVLGSPKRKQ HDVTQSKLFS LKTSMIQHDK LNLGPAIKNR YSSNIQMQAA NNLNQENSRP LQEKNTESSE NMTKFPSSRG KSTVSLNKNK QLSNSAVFVV SNNNLGVVKS AVDFQMKGKE NLTGKGITQI LGKSHGSLKL KRQPHIFESD GE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com