

Product Data Sheet

Recombinant Human METTL2A Protein

Catalog #	Source	Reactivity	Applications
CST1198	E. coli	H	
Description	Recombinant Human METTL2A Protein is produced by E. coli expression system and the target gene encoding Glu51-Ile170 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	METTL2A		
Alternative Names	METTL2; Methyltransferase-like protein 2A		
Entrez Gene	339175		
SwissProt	Q96IZ6		
Purity			
Chemical Structure	ERKVQENSIQ RVCQEKQVDY EINAHKYWND FYKIHENGFF KDRHWLFTEF PELAPSQNQN HLKDWFLLENK SEVPECRNNE DGPGLIMEEQ HKCSSKSLEH KTQTLPVEEN VTQKISDLEI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com