

Product Data Sheet

Recombinant Human TSGA10IP Protein

Catalog #	Source	Reactivity	Applications
CST1199	E. coli	H	
Description	Recombinant Human TSGA10IP Protein is produced by E. coli expression system and the target gene encoding Cys161-Thr310 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TSGA10IP		
Alternative Names	Testis-specific protein 10-interacting protein; Tsga10-interacting protein		
Entrez Gene	254187		
SwissProt	Q3SY00		
Purity			
Chemical Structure	CWKTEAQNLK ARQQLGAWGG VSIPTGKGEL GSEPPSGLQL PGRRPGSGSA SDKQVQLQSL GAEEAERGLS SGVLPQRPRR GSISEEEQFS EATEEAEEGE HRTPCRRRAG CQRKGQISGE EASDEGEVQG QSQGSSPSFN NLRRRQWRKT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com