

Product Data Sheet

Recombinant Human MTFR1L Protein

Catalog #	Source	Reactivity	Applications
CST1205	E. coli	H	
Description	Recombinant Human MTFR1L Protein is produced by E. coli expression system and the target gene encoding Glu191-Arg260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MTFR1L		
Alternative Names	FAM54B; Mitochondrial fission regulator 1-like		
Entrez Gene	56181		
SwissProt	Q9H019		
Purity			
Chemical Structure	EETEVEVPEL PSVPLLCAS PECKPEHKA ACSSEEDDC VLSKASSFA DMMGILKDFH RMKQSQDLNR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com