

Product Data Sheet

Recombinant Human EVPLL Protein

Catalog #	Source	Reactivity	Applications
CST1207	E. coli	H	
Description	Recombinant Human EVPLL Protein is produced by E. coli expression system and the target gene encoding Asp11-Gln260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EVPLL		
Alternative Names	Envoplakin-like protein		
Entrez Gene	645027		
SwissProt	A8MZ36		
Purity			
Chemical Structure	DILETQKRLQ QDRLNSEQSQ ALQHQQETGS SLKEAEVLLK DLFLDVDKAR RLKHPQAEET EKDIEQLHER VTQECAEYCA LYEKMVLPPR RGIQGR LGTR AGAETEAGLR RPVWAGHGGA GGTDRGAQHR AEGDQRPRRA AAEPGGAGCR HHPEPIPRPT EGGVVARAEP GQPVHALQGC TWQLSALAEQ QRRILQQDWS DLMADPAGVR REYEHFKQHE LLSQEQSVNQ LE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com