

Product Data Sheet

Recombinant Human WDR89 Protein

Catalog #	Source	Reactivity	Applications
CST1209	E. coli	H	
Description	Recombinant Human WDR89 Protein is produced by E. coli expression system and the target gene encoding Ala141-Lys270 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	WDR89		
Alternative Names	C14orf150; WD repeat-containing protein 89		
Entrez Gene	112840		
SwissProt	Q96FK6		
Purity			
Chemical Structure	ALLVFWDARM NSQNLSTTKD SLGAYSETHS DDVTQVRFP SNPNMVVSGS SDGLVNVFDI NIDNEEDALV TTCNSISSVS CIGWSGKGYK QIYCMTHDEG FYWWDLNHLD TDEPVTRLNI QDVREVVNMK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com