

Product Data Sheet

Recombinant Human HGH1 Protein

Catalog #	Source	Reactivity	Applications
CST1210	E. coli	H	
Description	Recombinant Human HGH1 Protein is produced by E. coli expression system and the target gene encoding Glu261-Pro380 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM203A		
Alternative Names	BRP16; C8orf30A; Protein FAM203A; Brain protein 16		
Entrez Gene	51236		
SwissProt	Q9BTY7		
Purity			
Chemical Structure	EEMERLPVDL QYLPPDKQRE PDADIRKMLV EAIMLLTATA PGRQQVRDQG AYLILRELHS WEPEPDVRTA CEKLIQVLIG DEPERGMENL LEVQVPEDVE QQLQQLDCRE QEQLERELAP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com