

Product Data Sheet

Recombinant Human LRRC75A Protein

Catalog #	Source	Reactivity	Applications
CST1211	E. coli	H	
Description	Recombinant Human LRRC75A Protein is produced by E. coli expression system and the target gene encoding Pro161-Ala340 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRRC75A		
Alternative Names	C17orf76; FAM211A; Leucine-rich repeat-containing protein 75A; Leucine-rich repeat-containing protein FAM211A		
Entrez Gene	388341		
SwissProt	Q8NAA5		
Purity			
Chemical Structure	PQACLKAVLA GSPPDNTVDL SGIPLTSRDL ERVTSYLQRC GEQVDSVELG FTGLTDDMVL QLLPALSTLP RLTTLALNGN RLTRAVLRDL TDILKDPSKF PNVTWIDLGN NVDIFSLPQP FLLSLRKRSP KQGHLP TILE LGEGPGSGEE VREGTVGQED PGGGPVAPAE DHHEGKETVA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com