

Product Data Sheet

Recombinant Human CCDC15 Protein

Catalog #	Source	Reactivity	Applications
CST1216	E. coli	H	
Description	Recombinant Human CCDC15 Protein is produced by E. coli expression system and the target gene encoding Gln101-Gln530 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CCDC15		
Alternative Names	Coiled-coil domain-containing protein 15		
Entrez Gene	80071		
SwissProt	Q0P6D6		
Purity			
Chemical Structure	QQLQKSYERA QKEGSIAMQS SATHLTSKRT SVFPNNLNVA IGSSRLPPSL MPGDGIEDEE NQNELFQQQA QALSETMKQA RHRLASFKTV IKKKGSVFPD DGRKSFLTRE EVLSRKPAST GINTGIRGEL PIKVHQGLLA AVPYQNYMEN QELDYEEP DY EESSSLVTDE KGKEDLFGRG QQDQQAIHSE DKNKPFSRVQ KVKFKNPLFV LMEEEEQKQL HFEGLDILP EA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com