

Recombinant Human GRASP65 Protein

Catalog #	Source	Reactivity	Applications
CST1218	E. coli	H	
Description	Recombinant Human GRASP65 Protein is produced by E. coli expression system and the target gene encoding Arg201-Leu430 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GORASP1		
Alternative Names	GOLPH5; GRASP65; Golgi reassembly-stacking protein 1; Golgi peripheral membrane protein p65; Golgi phosphoprotein 5; GOLPH5; Golgi reassembly-stacking protein of 65 kDa; GRASP65		
Entrez Gene	64689		
SwissProt	Q9BQQ3		
Purity			
Chemical Structure	RIPTQPPSYH KKPPGTTPPS ALPLGAPPPD ALPPGPTPED SPSLETGSRQ SDYMEALLQA PGSSMEDPLP GPGSPSHSAP DPDGLPHFME TPLQPPPPVQ RVMDPGFLDV SGISLLDNSN ASVWPSLPSS TELTTAVST SGPEDICSSS SSHERGGEAT WSGSEFEVSF LDSPGAQAQA DHLPLQLTPD SLTSAASPED GLSAELLEAQ AEEEPASTEG LDTGTEAEGL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com