

Recombinant Human CDK5RAP2 Protein

Catalog #	Source	Reactivity	Applications
CST1219	E. coli	H	
Description	Recombinant Human CDK5RAP2 Protein is produced by E. coli expression system and the target gene encoding Glu271-Lys480 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDK5RAP2		
Alternative Names	CEP215; KIAA1633; CDK5 regulatory subunit-associated protein 2; CDK5 activator-binding protein C48; Centrosome-associated protein 215		
Entrez Gene	55755		
SwissProt	Q96SN8		
Purity			
Chemical Structure	ERETAAQME HQKERNSFEE RIQALEEDLR EKEREIATEK KNSLKRDKAI QGLTMALKSK EKKVEELNSE IEKLSAAFAK AREALQKAQT QEFQGSSEDYE TALSGKEALS AALRSQNLTK STENHRLRRS IKKITQELSD LQQERERLEK DLEEAHREKS KGDCTIRDLR NEVEKLRNEV NEREKAMENR YKSLLSESNN KLHNQEQVIK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com