

Recombinant Human ZNF395 Protein

Catalog #	Source	Reactivity	Applications
CST1223	E. coli	H	
Description	Recombinant Human ZNF395 Protein is produced by E. coli expression system and the target gene encoding Lys11-His220 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZNF395		
Alternative Names	HDBP2; PBF; Zinc finger protein 395; HD-regulating factor 2; HDRF-2; Huntington disease gene regulatory region-binding protein 2; HD gene regulatory region-binding protein 2; HDBP-2; Papillomavirus regulatory factor 1; PRF-1; Papillomavirus-binding factor		
Entrez Gene	55893		
SwissProt	Q9H8N7		
Purity			
Chemical Structure	KRSLGARVL GPSASEGPSA APPSEPLLEG AAPQPFTTSD DTPCQEQPKE VLKAPSTSGL QQVAFQPGQK VYVWYGGQEC TGLVEQHSWM EGQVTWVLE QKLQVCCRVE EVWLAELQGP CPQAPPLEPG AQALAYRPVS RNIDVPKRKS DAVEMDEMMA AMVLTSLSGS PVVQSPPGTE ANFSASRAAC DPWKESGDIS DSGSSTTSGH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com