

Product Data Sheet

Recombinant Human RBM12B Protein

Catalog #	Source	Reactivity	Applications
CST1235	E. coli	H	
Description	Recombinant Human RBM12B Protein is produced by E. coli expression system and the target gene encoding Gln481-Glu680 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RBM12B		
Alternative Names	RNA-binding protein 12B; RNA-binding motif protein 12B		
Entrez Gene	389677		
SwissProt	Q8IXT5		
Purity			
Chemical Structure	QEFGVNF SVM SSEKMQARSQ SRERGDHSHL FDSKDPPIYS VGAFENFRHQ LEDLRQLDNF KHPQRDFRQP DRHPPEDFRH SSED FRFPPE DFRHSPEDFR RPRED FRRP SEED FRRPWE EDFRRPPEDD FRHPREEDWR RPLEEDWRRP LEED FRRSPT EDFRQLPEED FRQPPEEDLR WLPEED FRRP PEEDWRRPPE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com