

Product Data Sheet

Recombinant Human NBEAL1 Protein

Catalog #	Source	Reactivity	Applications
CST1250	E. coli	H	
Description	Recombinant Human NBEAL1 Protein is produced by E. coli expression system and the target gene encoding Glu1301-Ala1590 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NBEAL1		
Alternative Names	ALS2CR16; ALS2CR17; Neurobeachin-like protein 1; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 16 protein; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 17 protein		
Entrez Gene	65065		
SwissProt	Q6ZS30		
Purity			
Chemical Structure	EDTKKNSDEK TDEEKITSFA SANVSSDQWS LEDRHSLSN TPLFPEDSSV GELSFKSENQ EEFWHSNPSH LSLDLGIDS CEMSDSGSQV PDSLPSTPSP VESTKSFSVH SDRESSITND MGFSDDFSLL ESQERCEEEL LQLLTHILNY VMCKGLEKSD DDTWIERGQV FSALSKPGIS SELLRPSDEI KLTLQKMLE WAISENREAK TNPVTAENAF RLVLIQDFL QS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com