

Product Data Sheet

Recombinant Human GPATCH8 Protein

Catalog #	Source	Reactivity	Applications
CST1258	E. coli	H	
Description	Recombinant Human GPATCH8 Protein is produced by E. coli expression system and the target gene encoding Asp11-Gly240 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GPATCH8		
Alternative Names	GPATC8; KIAA0553; G patch domain-containing protein 8		
Entrez Gene	23131		
SwissProt	Q9UKJ3		
Purity			
Chemical Structure	DRDFQGNHFD QYEEGHLEIE QASLDKPIES DNIGHRLLQK HGWKLGQGLG KSLQGRTDPI PIVVKYDVMG MGRMEMELDY AEDATERRRV LEVEKEDTEE LRQKYKDYVD KEKAIKALE DLRANFYCEL CDKQYQKHQE FDNHINSYDH AHKQRLKDLK QREFARNVSS RSRKDEKKQE KALRRLHELA EQRKQAEAP GSGPMFKPTT VAVDEEGGED DKDESATNSG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com