

Recombinant Human DIS3 Protein

Catalog #	Source	Reactivity	Applications
CST1260	E. coli	H	
Description	Recombinant Human DIS3 Protein is produced by E. coli expression system and the target gene encoding Ala191-Leu330 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DIS3		
Alternative Names	KIAA1008; RRP44; Exosome complex exonuclease RRP44; Protein DIS3 homolog; Ribosomal RNA-processing protein 44		
Entrez Gene	22894		
SwissProt	Q9Y2L1		
Purity			
Chemical Structure	AFTCEEYVKS LTANPELIDR LACLSEEGNE IESGKIIFSE HLPLSKLQQG IKSGTYLQGT FRASRENYLE ATVWIHGDNE ENKEILQGL KHLNRAVHED IVAVELLPKS QWVAPSSVVL HDEGQNEEDV EKEEETERML		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com