

Recombinant Human HECW1 Protein

Catalog #	Source	Reactivity	Applications
CST1262	E. coli	H	
Description	Recombinant Human HECW1 Protein is produced by E. coli expression system and the target gene encoding Val61-Cys550 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HECW1		
Alternative Names	KIAA0322; NEDL1; E3 ubiquitin-protein ligase HECW1; HECT, C2 and WW domain-containing protein 1; NEDD4-like E3 ubiquitin-protein ligase 1; hNEDL1		
Entrez Gene	23072		
SwissProt	Q76N89		
Purity			
Chemical Structure	VTIPRSTSDT DLVTSDSRST LMVSSSYYSI GHSQDLVIHW DIKEEVDAGD WIGMYLIDEV LSENFLDYKN RGVNGSHRGQ IWKIDASSY FVEPETKICF KYYHGVSGAL RATTPSVTVK NSAAPIFKSI GADETVQGQG SRRLISFSL SDFQAMGLKKG MFFNPDPYLK ISIQPGKHSI FPALPHHGQE RRSKIIGNTV NPIWQAEQFS FVSLPTDVLE IEVKDKFAKS RP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com