

Recombinant Human CFAP45 Protein

Catalog #	Source	Reactivity	Applications
CST1270	E. coli	H	
Description	Recombinant Human CFAP45 Protein is produced by E. coli expression system and the target gene encoding His61-Glu460 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CCDC19		
Alternative Names	NESG1; Coiled-coil domain-containing protein 19, mitochondrial; Nasopharyngeal epithelium-specific protein 1		
Entrez Gene	25790		
SwissProt	Q9UL16		
Purity			
Chemical Structure	HTLQKLTAL GLDRKPETIQ LITRDMVREL IVPTEDPSGE SLIISPEEFE RIKWASHVLT REELEARDQA FKKEKEATMD AVMTRKKIMK QKEMVWNNNK KLSDLEEVAK ERAQNLLQRA NKLRMEQEEE LKDMSKIILN AKCHAIRDAQ ILEKQQIQKE LDTEEKRLDQ MMEVERQKSI QRQEELERKR REERIRGRRQ IVEQMEKNQE ERSLLAEQRE QEKEQMLEYM EQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com