

Product Data Sheet

Recombinant Human CEP170B Protein

Catalog #	Source	Reactivity	Applications
CST1271	E. coli	H	
Description	Recombinant Human CEP170B Protein is produced by E. coli expression system and the target gene encoding Met281-Ala450 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CEP170B		
Alternative Names	FAM68C; KIAA0284; Centrosomal protein of 170 kDa protein B; Centrosomal protein 170B; Cep170B		
Entrez Gene	283638		
SwissProt	Q9Y4F5		
Purity			
Chemical Structure	MKIKDHITKF SLRQRRPPGK EATPGEMVSA ETKVADWLQVQ NDPSSLHRVG PGDDRHSKTS DLPVHTRTLK GHKHEDGTQS DSEDPLAKAA SAAGVPLEAS GEQVRLQRQI KRDPQELLHN QQAFVIEFFD EDTPRKKRSQ SFTHSPSGDP KADKRRGPTP ADRDRPSVPA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com