

Product Data Sheet

Recombinant Human BEND2 protein

Catalog #	Source	Reactivity	Applications
CST1273	E. coli	H	
Description	Recombinant Human BEND2 protein is produced by E. coli expression system and the target gene encoding Ile641-Gly790 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	BEND2		
Alternative Names	CXorf20; BEN domain-containing protein 2		
Entrez Gene	139105		
SwissProt	Q8NDZ0		
Purity			
Chemical Structure	IPEGMREPST DNPEEPGEAW SYFGRPWRNI RMPCSVLTLA KTKSCASLSA RYLIQKLFTK DVLVQSNVYG NLKHGLCALD PNKISALREF LQENYPICDL SENG RDWKSC VTSINSGIRS LRHDVRRAEA RSQSLPAVTP PELEQESKPG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com