

Product Data Sheet

Recombinant Human ARAP3 Protein

Catalog #	Source	Reactivity	Applications
CST1275	E. coli	H	
Description	Recombinant Human ARAP3 Protein is produced by E. coli expression system and the target gene encoding Leu71-Glu240 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ARAP3		
Alternative Names	CENTD3; Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3; Centaurin-delta-3; Cnt-d3		
Entrez Gene	64411		
SwissProt	Q8WWN8		
Purity			
Chemical Structure	LDPKSDSAME PSPSPAPQAA PPKVPKPRV VFGGLSGPAT TQRPGLSPAL GGPGVSRSP PSPRPPLPT SSSEQSSALN TVEMMPNSIY FGLDSRGRAQ AAQDKAPDSS QISAPTPALR PTTGTVHIMD PGCLYYGVQP VGTPGAPDRR ESRGVCQGRA EHRLSRQDLE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com