

Product Data Sheet

Recombinant Human ASB15 Protein

Catalog #	Source	Reactivity	Applications
CST1276	E. coli	H	
Description	Recombinant Human ASB15 Protein is produced by E. coli expression system and the target gene encoding His11-Val430 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ASB15		
Alternative Names	Ankyrin repeat and SOCS box protein 15; ASB-15		
Entrez Gene	142685		
SwissProt	Q8W XK1		
Purity			
Chemical Structure	HLTSYDIQLS IQESIEASKT ALCPERFVPL SAQNRKLVEA IKQGHIPQLQ EYVKYKYAMD EADEKGWFPL HEAVVQPIQQ ILEIVLDASY KTLWEFKTCD GETPLTLAVK AGLVENVRTL LEKGVWPNTK NDKGETPLLI AVKKGSYDMV STLIKHNTSL DQPCVKRWSA MHEAAKQGRK DIVALLLKHG GNVHLRDGFG VTPLGVAAEY GHCDVLEHLI HKGGDVLALA DD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com