

Product Data Sheet

Recombinant Human RNF115 Protein

Catalog #	Source	Reactivity	Applications
CST1280	E. coli	H	
Description	Recombinant Human RNF115 Protein is produced by E. coli expression system and the target gene encoding Gln221-Asn290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RNF115		
Alternative Names	ZNF364; E3 ubiquitin-protein ligase RNF115; RING finger protein 115; Rabring 7; Zinc finger protein 364		
Entrez Gene	27246		
SwissProt	Q9Y4L5		
Purity			
Chemical Structure	QVDMGLECPV CKEDYTVEEE VRQLPCNHFF HSSCIVPWLE LHDTCPVCRK SLNGEDSTRQ SQSTEASASN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com