

Product Data Sheet

Recombinant Human UBXL6 Protein

Catalog #	Source	Reactivity	Applications
CST1291	E. coli	H	
Description	Recombinant Human UBXL6 Protein is produced by E. coli expression system and the target gene encoding Ala31-Glu110 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UBXL6		
Alternative Names	UBXD1; UBXD2; UB domain-containing protein 6; UB domain-containing protein 1		
Entrez Gene	80700		
SwissProt	Q9BZV1		
Purity			
Chemical Structure	AHKEKPNQPA PRPPRQGPTN EAQMAAAAAL ARLEKQKSRA WGPTSQDTIR NQVRKELQAE ATVSGSPEAP GTNVVSEPRE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com