

Product Data Sheet

Recombinant Human RBM42 Protein

Catalog #	Source	Reactivity	Applications
CST1292	E. coli	H	
Description	Recombinant Human RBM42 Protein is produced by E. coli expression system and the target gene encoding Glu391-Lys470 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RBM42		
Alternative Names	RNA-binding protein 42; RNA-binding motif protein 42		
Entrez Gene	79171		
SwissProt	Q9BTD8		
Purity			
Chemical Structure	EVNDDILARA FSRFPSFLKA KVIRDKRTGK TKGYGFSVSK DPSDYVRAMR EMNGKYVGSR PIKLRKSMWK DRNLDVVRKK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com