

Recombinant Human TRIM41 Protein

Catalog #	Source	Reactivity	Applications
CST1297	E. coli	H	
Description	Recombinant Human TRIM41 Protein is produced by E. coli expression system and the target gene encoding Glu71-Glu390 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRIM41		
Alternative Names	RINCK; E3 ubiquitin-protein ligase TRIM41; RING finger-interacting protein with C kinase; RINCK; Tripartite motif-containing protein 41		
Entrez Gene	90933		
SwissProt	Q8WV44		
Purity			
Chemical Structure	EEVEAVGAGA GWDTPMRDED YEGDMEEVE EEEGVFWTS GMSRSSWDNM DYVWEEEDDEE EDLDYYLGDM EEEDLRGEDE EDEEEVLEEV EEEDLDPVTP LPPPPAPRRRC FTCPQCRKSF PRRSFRPNLQ LANMVQVIRQ MHPTPGRGSR VTDQGICPKH QEALKLFCEV DEEAICVVCR ESRSHKQHSV VPLEEVVQEY KAKLQGHVEP LRKHLEAVQK MKAKEERRVT EL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com