

Product Data Sheet

Recombinant Human AFAP1L2 Protein

Catalog #	Source	Reactivity	Applications
CST1302	E. coli	H	
Description	Recombinant Human AFAP1L2 Protein is produced by E. coli expression system and the target gene encoding Pro641-Val740 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	AFAP1L2		
Alternative Names	KIAA1914; XB130; Actin filament-associated protein 1-like 2; AFAP1-like protein 2		
Entrez Gene	84632		
SwissProt	Q8N4X5		
Purity			
Chemical Structure	PPVKDRLRVT SAEIKLGKNR TEAEVKRYTE EKERLEKKKE EIRGHLAQLR KEKRELKETL LKCTDKEVLA SLEQKLKEID EECRGEESRR VDLELSIMEV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com