

Product Data Sheet

Recombinant Human NELFCD Protein

Catalog #	Source	Reactivity	Applications
CST1306	E. coli	H	
Description	Recombinant Human NELFCD Protein is produced by E. coli expression system and the target gene encoding Asp11-Leu220 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NELFCD		
Alternative Names	NELFD; TH1; TH1L; Negative elongation factor C/D; NELF-C/D; TH1-like protein		
Entrez Gene	51497		
SwissProt	Q8IXH7		
Purity			
Chemical Structure	DEDYYGSAAE WGDEADGGQQ EDDSGEGEDD AEVQQECLHK FSTRDYIMEP SIFNTLKRYF QAGGSPENVI QLLSENYTAV AQTVNLLAEW LIQTGVEPVQ VQETVENHLK SLLIKHFDPR KADSIFTEEG ETPAWLEQMI AHTTWRDLFY KLAEAHDPDCL MLNFTVKLIS DAGYQGEITS VSTACQQLEV FSRVLRTSLA TILDGGEENL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com