

Product Data Sheet

Recombinant Human WWC1 Protein

Catalog #	Source	Reactivity	Applications
CST1307	E. coli	H	
Description	Recombinant Human WWC1 Protein is produced by E. coli expression system and the target gene encoding Thr811-Arg1050 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	WWC1		
Alternative Names	KIAA0869; Protein KIBRA; HBeAg-binding protein 3; Kidney and brain protein; KIBRA; WW domain-containing protein 1		
Entrez Gene	23286		
SwissProt	Q8IX03		
Purity			
Chemical Structure	TDAVSALLEQ TAVELEKRQE GRSSTQTLED SWRYEETSEN EAVAEEEEEE VEEEEGEEDV FTEKASPDMD GYPALKVDKE TNTETPAPSP TVVRPKDRRV GTPSQGPFLR GSTIIRSKTF SPGPQSQYVC RLNRSDSDSS TLSKKPPFVR NSLERRSVRM KRPSSVKSLR SERLIRTSLD LELDLQATRT WHSQLTQEIS VLKELKEQLE QAKSHGEKEL PQWLREDERF RL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com