

Product Data Sheet

Recombinant Human SPATA24 Protein

Catalog #	Source	Reactivity	Applications
CST1308	E. coli	H	
Description	Recombinant Human SPATA24 Protein is produced by E. coli expression system and the target gene encoding Glu51-Gln160 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SPATA24		
Alternative Names	Spermatogenesis-associated protein 24; Testis protein T6441 homolog		
Entrez Gene	202051		
SwissProt	Q86W54		
Purity			
Chemical Structure	EEFQAVEKKL VEEKAAHAKT KVLLAKEEEK LQFALGEVEV LSKQLEKEKL AFEKALSSVK SKVLQESSKK DQLITKCNEI ESHIKQEDI LNGKENEIKE LQQVISQQKQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com